

Product datasheet for **RR213934L3V**

Srgn (NM_020074) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Srgn (NM_020074) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Srgn
Synonyms:	Pgsg
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020074
ORF Size:	537 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213934).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020074.2 , NP_064459.1
RefSeq Size:	1320 bp
RefSeq ORF:	540 bp
Locus ID:	56782
UniProt ID:	P04917
Cytogenetics:	20q11
Gene Summary:	precursor for the core protein of the yolk sac carcinoma proteoglycan, PG19, component of extracellular matrices, cell membranes, and some intracellular organelles [RGD, Feb 2006]


[View online »](#)