

## Product datasheet for **RR213941L4V**

### Elob (NM\_031129) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Elob (NM_031129) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Elob
Synonyms:	Tceb2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031129
ORF Size:	354 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213941).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031129.3</a> , <a href="#">NP_112391.1</a>
RefSeq Size:	484 bp
RefSeq ORF:	357 bp
Locus ID:	81807
UniProt ID:	<a href="#">P62870</a>
Cytogenetics:	10q12
Gene Summary:	positive regulatory subunit of the general transcription factor SIII heterotrimer; acts to increase the rate of RNA polymerase II transcription; regulates RNA polymerase II elongation complex activity [RGD, Feb 2006]



[View online »](#)