

Product datasheet for **RR213961**

LOC499796 (NM_001047951) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LOC499796 (NM_001047951) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LOC499796
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

```
>RR213961 representing NM_001047951
Red=Cloning site Blue=ORF Green=Tags(s)
```

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCTGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**

ATGATAGGCGAGTACAGAAAGCTTCAAGGGCAGACACCCACCTTCTCCAGCAGTCATTACTGATTCTGTACCTG
 GTTTGAGCAATCCAGCCGAGTGCCTCCAAAGCCTCCTGGAGCCTACACCATTGTCTATGTACTTGTCTGCT
 TTCTGCCCTAAGGCCTCTGCAGTATAGGGGGTGGCTAGGCGCAGAGGTAATGTCTTTGATCTCTGCCCC
 TTGCCACCTCCGGGAATCGAGAAAGCTGACTACAGGGTCTCTGAGCTACCAAGTCAAAGGTATCAGGTACCA
 GGGAAAGGTATTGCTACCAAGTTGCCAGCTGCATCTTTTGGGAGAGCGGGCCCTGCATCTCACCTGGACAG
 CACAGCAGAGCTAACCTTGGTAGAGTGGTGCAGGTGAGCCAGCCCGGAGGGCAAAAAAGTGGTAGAGC
 GGCCCCACCATTCATCAGCGCTACCATCTACAACCAAAATAGTTCTCAACCCGAGTGAGGCTTCATGAAC
 TCACCCCAACCTGCCGAATGTTTGTCTGGTCAGTCATCCACAGCTGCTGTGACTTCTGGGTGCAGTT
 TTTTCAACACACACATTGCCTTAATTTTAGTGTGCATGTACAGTTTCATCATCTAGTGATTCTAAGCAG
 AGTGGAGTGCACAAGGAGCCAGGAAGCGGGGCAGCTGCACCTAAAAACAGTGACACAGAACTCCT
 CCAGATTGAAGCAGACAGCATTAATCAGGAAGTCCCTTGAGGTGAGATACCTGGCATCAGGAAGACTTATA
 CTATGAGAGAAAGCATTCAGGATACCTGAGCTGCTTGGGCAATACTTACGGGATTCTTCCAGAAGCTCTA
 AAATCCGACTTAACTGAGTACTGTGCACTTGTCAAAGCTGTGGATGAGAACAATGACCTGTGCCAGTGTG
 GCATATTTGGGGCTGATTTCTGCTCTCTTTCAGACTTTTACCTTACACCAGCAGGGACTGAGTCTCTC
 TCTACCTCACAGAACAAACAAAGGTGCTGTGTTGCTGCCAATTTTACCACAAAGCATCCAGGAGGTCCA
 GAATTGAGAAGAAGAACCCTCTTTCATTTTACCTTCATCAAGGTGGAAAGTGAAGGAAAAGTGGTGCAGT
 TGGACACCGAGAGTCTCGGAGGGTGTTCTGAAGATAACACACCACCAGGGAGACAGTTTGGCTGTACATG
 CACAAGGAGAGCCTCATCATGGACAGCAGTGGGTGAACTGAGACAGACTGCTCTCAAAGAAGTCAGTGT
 GTGAAAGAGAAGGATGGAAGTGAAGAAATCCTCCAGGGGCAGAGGCTGGATAACGGCCATCTTTAATAGGA
 GCGGGAACGCAGGAAAAAGAAAGAACATCTAGGAAGGACACCATTGAGATATCCGTCAAATAGCAAACAGGA
 AAAAGAAGGGGGAGGAGGAAGAGGAGGAGAAGGAGAAGAGAAGAGGAGAAGAGAAGGCAAGAGGAG
 AAGAGGAGGAGGAAAGGAAAGAAAGAAAGAAAGAAAGAAAGAAAGAAAGAAAGAAAGAAAGAAAGAA
 AGAAAGAAAGAAAGGAAGAAGGAAAGGGAAGGAAGGAAGGAAGGAAGGAAGGAAGGAAGGAAGGAAGAA
 GAAAGAAAGAAAGAAAGGGAAGGAAGGAAGGAAGGAAGGAAGGAAGGAAGGAAGGAAGGAAGGAAGAA
 CCAAGAAACAAAAACACAGAGAGAAAGACCCGGGTCTCCTGCCCTGTATGCCGGTCAAGTTGCAGCTGAGT
 CTGTGAAGTCTGTGAAGGTCTTTGTGTTTGGAGTTGCAGGAAATCCTCTTGATTTACGGAAAGTTCTCT
 CTACGCTGGACTGAGCAGCTGGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA

Protein Sequence:


>RR213961 representing NM_001047951
Red=Cloning site Green=Tags(s)

MKAVQKLQGQNPSPAVITDSVPGLSNPAECLQSLEPTPLSMYLLL SALRPLQYRGWPRRRGNVFDLCP
 LPPPGIEKADYRVL SYQSVSGTRERYCYQLPAASFGRAGPASHLDSTAE LTLVEVVQVSQPGGQKRGRA
 GPTIHQRYHLQPNSSSTRVRLHELTPNPAGMFAQSSTA AVTSWVQFFSTPHIALILVCMYSSSSSDSKQ
 SGVRNKEPGSGARRPKNTVTQNSSRVSTEHNQEVPLRSDTWHQEDLYYERKHSGYLTC LGNTYRILPERL
 KSDLTEYCALVKAVDENNDLCQCGIFGADIPALFQTFTLHQQGLSLSLPHRTNKGAVLLPIFTTKHPGGP
 ELRRRTLHFHTFIKVESEGKLVQLDTESLGGCSEDNTPPGRQFGCTCTRASSWTAVGELRQTALKEVSV
 VKEKDGLKNPPRGRGWITAI FNRSNAGKEEHLGRTPLYPSNSKQEKEGGGGRGGEGEREERGEEKAKE
 KRRRKGGKKERKKERKKERKKERKEEGKRKEGRKEGRKKERKKERKEGRKEGRKEGRKKRRN
 KNKTNREKRPGPSALYAGOVAESVKS VKVFVFGVAGNPLGFTESSLRWTEQLV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001047951

ORF Size: 1845 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001047951.1](#), [NP_001041416.1](#)

RefSeq Size: 1848 bp


RefSeq ORF: 1848 bp

Locus ID: 499796

Cytogenetics: 3q12

MW: 68.9 kDa

Product images:



Circular map for RR213961