

Product datasheet for **RR213967**

Ctsql2 (NM_001002813) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ctsql2 (NM_001002813) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ctsql2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR213967 representing NM_001002813 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGACTGCTGCTCTTCTCTGATCATCCTGTGCTGGGAGTTGTGTCAAGTGCTTCAGCATTTAATCTCA
 GTTTGGATGTCCAATGGCAAGAAATGGAAGATGAAATATGAAAAATTATACAGTCCGGAGGAAGAACTCCT
 GAAGAGAGTGGTATGGGAAGAAAATGTGAAAAAGATTGAAGTGCATAACAGGGAGAATTCCTTGGGAAG
 AATACCTATATTATGAAATAAAATAATTTTGTGACCTGACTGATGAAGAATTCAGGACATGATAACTG
 GTATTACGTTGCCAATTAACAACACAATGAAAAGTCTCTGGAAACGTGCACTTGGTAGTCCTTTTCCTAA
 TTCTTGGTATTGGCGAGATGCTTTGCCCAAATCCATTGATTGGAGAAAGGAAGGCTATGTGACTCGTGTG
 AGGGAACAGGGAAGTGAAGTCTTGTGGGCTTTTCTGTAGCTGGTGCATAGAAGGACAAATGTTCA
 AGAAAACAGGCAAACTGACCCCACTGAGTGTCCAGAATCTTGTGGACTGTTCTAAACCTCAAGGCAATAA
 AGGCTGTCGTGGAGGTACTACATACAATGCATTCCAGTATGTTTTGCAAAATGGAGGTTTGGAGTCTGAG
 GCAACCTATCCATATAAAGGAAAAGAAGGACTATGCAAGTACAATCCTAAGAATGCATATGCTAAAATCA
 CAAGATTTGTGCTCTCCAGAAGATGAAGATGCCTAATGGATGCTTTAGCAACTAAAGGACCTGTTGC
 TGCTGGAATTCATGTTGTCTACAGTTCTTTAAGGTTCTACAAAAAAGGTATATATCATGACCCAAAGTGC
 AACATCGTGTGAATCATGCAGTTCTCGTGGTTGGCTATGGATTTGAGGGAAATGAAACAGATGGCAATA
 ACTACTGGCTGATCAAGAACAGCTGGGGTAAACAGTGGGGATTGAAGGCTATATGAAGATCGCCAAAGA
 CAGGAACAACCACTGTGGAATTGCCACATTTGCCCAATACCCTATTGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RR213967 representing NM_001002813

Red=Cloning site Green=Tags(s)



MTAALFLIILCLGVVSGASAFNLSLDVQWQEWKMYEKLYSPEEELLKRVVWEENVKKIELHNRENSLGK
NTYIMEINNFDALTDEEFKDMITGITLPINNTMKS LWKRALGSPFPNSWYWRDALPKSIDWRKEGYVTRV
REQQCKCKSCWAFPVAGIAEQMFKKTGKLTPLSVQNLVDCSKPGQNKGRGGTIVYNAFYQLVQNGGLESE
ATPYQKGEKGLCKYNPKNAYAKITRFVALPEDEVDLMDALATKPGVAAIGHVYSSLYFRYKKGIIYHEPKC
NNRVNHAVLVYGYGFEENETDGNMYLTKNSGWQOWGLGYMKIAKDRNNHCGIATFAOYPTV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM 001002813

ORF Size: 1029 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001002813.2, NP_001002813.2

RefSeq Size: 1319 bp


RefSeq ORF: 1032 bp

Locus ID: 408201

Cytogenetics: 17p14

MW: 39 kDa

Product images:



Circular map for RR213967