

Product datasheet for **RR213981L3V**

Hgh1 (NM_001007707) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hgh1 (NM_001007707) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hgh1
Synonyms:	Brp16; Fam203a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001007707
ORF Size:	1179 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213981).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001007707.1 , NP_001007708.1
RefSeq Size:	2288 bp
RefSeq ORF:	1182 bp
Locus ID:	315094
UniProt ID:	Q6AY79
Cytogenetics:	7q34


[View online »](#)