

Product datasheet for **RR213984L3V**

Abhd17c (NM_001100736) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Abhd17c (NM_001100736) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Abhd17c
Synonyms:	Fam108c1; RGD1308210
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001100736
ORF Size:	960 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213984).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001100736.1 , NP_001094206.1
RefSeq Size:	2208 bp
RefSeq ORF:	963 bp
Locus ID:	361601
UniProt ID:	B5DFK7
Cytogenetics:	1q31
Gene Summary:	Hydrolyzes fatty acids from S-acylated cysteine residues in proteins. Has depalmitoylating activity towards DLG4/PSD95.[UniProtKB/Swiss-Prot Function]



[View online »](#)