

Product datasheet for **RR213986L4V**

Defb4 (NM_022544) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Defb4 (NM_022544) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Defb4
Synonyms:	Defb2; Defb3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022544
ORF Size:	189 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213986).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022544.2 , NP_071989.1
RefSeq Size:	421 bp
RefSeq ORF:	192 bp
Locus ID:	64389
UniProt ID:	O88514
Cytogenetics:	16q12.5
Gene Summary:	cationic peptides with broad-spectrum antimicrobial activity; may mediate mucosal defenses of several organs [RGD, Feb 2006]


[View online »](#)