

Product datasheet for **RR213987**

Prrg2 (NM_001108482) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Prrg2 (NM_001108482) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Prrg2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR213987 representing NM_001108482
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGGAGCCATCCTTCCTGCTCTTAGTCTACATGGGATTGGCCACCTGCCTGGACACCTACCTCATA
 GAGAGCAAAGCCAAGTTCTTGACATCTTCCTGGATGCCCCAGAGGCCAGAGTTTCTTAGTTGGCCGCAG
 GCGTTTTCTCGTGCTAACCACTGGGACCTGGAGTTGCTGACCCCTGGGAACCTGGAACGGGAGTGTCTG
 GAAGAGAGGTGCTCTGGGAGGAGGCCGGGAGTACTTTGAGGACAACACACTGACGGAGCGTTTCTGGG
 AAAGCTATGCATACAACGGCAAGGAGGGCGTGGACGGGTGGATGTTGCAGGCCTGGCAGTGGGACTGAC
 CTCTGGGATCCTGCTCATCGTGTTGGCTGGCTTGGGAGCCTTTTGGTATCTGCATTATCGACGGCGCCGC
 CTCAGAGTCCAGGAGTCTGTCTGCAAGAGACCGGGCTCATCATCCCTTAAGTCCTCCGACGCCTCTGT
 CTCCACCGCCTCCCCAGGCCTTCCACCTATGAGCAGGCTTTAGCAGCCTCCGGGTGCACGACGCCCC
 TCCACCCCTTACTCCAGG

A**CGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR213987 representing NM_001108482
 Red=Cloning site Green=Tags(s)

MRSHPSLLL VYMGLATCLDTLPHREQSQVLDIFLDAPEAQSFVGRRRFPRANHWDELLTPGNLERECL
 EERCSWEEAREYFEDNTLTERFWESYAYNGKGRGRVDVAGLAVGLTSGILLIVLAGLGAFWYLHYRRRR
 LRVQESCLQETGLIIPSPPTPLSPPPPLPTYEQALAAAGVHDAPPPYSR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_001108482

ORF Size: 579 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001108482.1](#), [NP_001101952.1](#)

RefSeq Size: 779 bp


RefSeq ORF: 582 bp

Locus ID: 361570

Cytogenetics: 1q22

MW: 21.7 kDa

Product images:



Circular map for RR213987