

Product datasheet for **RR214018L4V**

Taf6 (NM_001044225) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Taf6 (NM_001044225) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Taf6
Synonyms:	TAF(II)80; TAFII70
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001044225
ORF Size:	1941 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214018).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001044225.1 , NP_001037690.1
RefSeq Size:	2127 bp
RefSeq ORF:	1944 bp
Locus ID:	288533
Cytogenetics:	12q11
Gene Summary:	TAFs are components of the transcription factor IID (TFIID) complex, PCAF histone acetylase complex and TBP-free TAFII complex (TFTC). TFIID is multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)