

Product datasheet for **RR214036**

Lrrc36 (NM_001004088) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lrrc36 (NM_001004088) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lrrc36
Synonyms:	Xlhsrf2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)


ORF Nucleotide
Sequence:

>RR214036 representing NM_001004088
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCGGAGCAATGGGATCTGGACGAGGAGTGCCTCCGCCGCTGGGGCGCTGACCCTCGAGCAGCCGG
 AACTGGTGGAGTCTCTCATTGCAGGGATCTTATGCTGGAAAAATCCACTCTATTGGTGTGCCTTCAG
 AAATTTCAAAAATCTCCGATCCTTAGATTTATCAAGGAATTTGATCACTAGTCTGAAGGGCATCCAGTAT
 CTATGTTCACTCCAAGATCTGAACCTATATTATAACAACATTCCTTCATTAGTGGAGGTGTCCCGACTAC
 AGCCCTTACCCTTCTCAAAGAAGTATTTGAGACTCAATCCTGTTGTAAGGAAAGATACAGATTATAG
 GCTCTTTGCTGTGTACACACTGCAATGCTGGAGAACTGGATGACCGAGCTGTCCGTGACAGTGAGAGG
 AGAGCGGCCAAGCTACATTTTCAGTCAGTGGGCAGCAGTAAAAATTTCTCCTAGAGGTGGAAAAAGCT
 CTAGGGAGAAGACAATGAAAACTGTGTGACAGATGAGGGCCCTGCATCACACGTCAGTTCTGAGGTTGA
 CACCAGAATGGAACAGACTCAAACAAGGACTTTTTATTTCCTTCCCAACCGGAAATAAAGGATTCC
 CTAGCAAGCATTGTGCAACCCAGGGCAGTGGCACACCAGCTCAGAAATTAGACGTTTTTCCACTGGAGA
 CACAGATGCAGGAAGCAACAAGAAGGGAGACAAGTGACCTCCACCAGGAAGATGGATTGGAACCTAGGCA
 TTAATCGCCTCATCAGGGCACAGTCCGATCTCCAGAGAAGATGACAAGAGATGGATACCGAGTATCTTTT
 TTAGACTAAGTCTTACAGTCTTCTCCAGTAAAGGACTTGATACCAAAACCTGATGCCTATCCGTGTA
 CCCATGATGCCTCTTGGGTAACGCTGGATGTGGGTGATTCTAACCAGATCCTTCCCTGCCAATTACC
 TTCTGGTGTGTCTGGACCATATGATAATCATTATTCTCAAACCTTTGCCTGCACGGAAGTCTTGTG
 AAAAGCCCTCAAAAAGGCAAGTACTATCGAGAATACAGCACAAGCCCTTCCAATGATAAAAAGCCGCTT
 CATCCCTTCTGTGGAGACTGCTAACATCCCTGTCAAACCTGACTGCAGCACAGGAAGGCCATTGAA
 GCTCAGTTCAGATCTCTATGCCACAACCCGTTTCAACAGTGACCATGCTCTGCTGGCCAATGTAGAGCAA
 CAGCTGTCCAGCCTCGGTGATTTTGACACAGCCTGGTTCCTTCCCAAGCAACCTGCCCTAGGAAATT
 CTCTGCGGACACTGCTGTGCCTCCTGGGACTCCAGAAAACAGAGAGATTCCAACCTAAGAGGTCTTAAG
 CCCTTCAAGGAGAGGATTCAAGCGGAAGGACAATATCCTTGCCAACCTGAATCTGAAACGTGGTTTCCAA
 GATGCTACAGGCAGCGAGGTGCAACCTTCTCTAGTGACCTGGCAGTTTGCATGGTTTATCAGGAAACC
 GTAGTCCCGGATCTCTGCCAGAACCCACGTAAGCCTGTCTCAGACAGCTCCTGGAGCTGGTGGAA
 TAAGCACTGGAATGGTTCTGGCTCCCTTCTCTCAACAAGAAGTTTCTTGGTCTGCCAGAGATTGCTC
 CTGTCTTGGTGGTGCCTGCACCTCCCAAGCAGTGGCATCGCTCGCTGCTTGACGACAAGACAGGGAGAG
 CCTTGCCTGGAGGGAGACGGAGCTAAAGGAAGCAGGGCAGCTGGTTCCCAAGTGTGGAAGTTTGAA
 GCAAAAACCTGGTCAAAGTGTGGAAGAAAACCTTATTTTGTGCAAAAAATTCAGCAGTTGGAGGGAGCT
 GCTGCCACCTCGGTTGTGAGTGGGCACCATCTCACACACAGATGATCTTCTGCACAAAAACCAACAGC
 TGACCATACAGGTGGCTTGCCTGAACAGGAGCTGACTCAGCTGAAAAGGCTGGAGGAGACAGTTGCTCT
 TCTCCATGAGAGTCAGAGTCCCTGGTGGTAACTAATGAGTATCTGCTGCAGCAGCTGAATAAAGAGCAG
 AAAGGTTATCCGGGAAATCGCTCCTGCCTCCTGAGAAGACTCATCCTTTGGGGAGATTGTCGCCCTTTG
 GGAAAAGCACTCTGTCTTCTCCTCACCAACGGTACACAACACGGGTAGTACCTAATACAGAGTGTCTC
 AGAAGCTGACCCTGAGCCTAGTCTGTGGAGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Plasmid Map:


ACCN: NM_001004088

ORF Size: 2271 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001004088.3](#), [NP_001004088.3](#)

RefSeq Size: 2410 bp

RefSeq ORF: 2274 bp

Locus ID: 361394

Cytogenetics: 19q12
MW: 84.1 kDa