

Product datasheet for **RR214039**

Gimap7 (NM_001024328) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gimap7 (NM_001024328) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gimap7
Synonyms:	Ian3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR214039 representing NM_001024328 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGCTGAACTGGATGATGACTCTCTGAGGATTGTTCTGGTAGGGAAAACAGGAAGTGAAAGAGCGCCA
CAGCAAACACCATCCTGGGGCAAAGATATTTACTTCTAGAATCGCACCCACGCTGTCACCCAGACCTG
TCAGAAGGCATCCCGCAGGTGGAAGGAAAGAGACCTCCTGGTTGTTGACACCCAGGACTCTTCGACACC
AAGGTGAACCTGGAACCAACAGCATTGAAATCAGCCGGTGTGTCCTTCAGTCCTGCCCCGGCCTCATG
CCATCATTCTTGTCTGCAACTGAATCGCTTCACGATAGAAGAGCAAGAGACCGTTACTAGGATCAAGGC
TATCTTTGGGAAGGCAGTCATGAAGTACCTGATCATCTTGTTACCCGCAAAGACGAGTTGGAGGACCAG
AACCTAAACGATTTTCATTGAAGACTCAGATACTAACCTTAAAAGCATCATCAAGGAGTGTGACAGCCGCT
ATCTGGCCATCAATAACAAAGCGGAGGGGGCTGAGGGGAAATGCAGGTGCAGGAAGTATGGGGTTTGT
AGAATCCCTGGTGCAGCAATGGCGGGCTCTACTTCTCTGATCCCATATACAAGTATGCAGAGCAAAGG
CTGAAGAAACAAGTGGGATCTTGAGGAAATTTACACTGACGTACTAGAGAAGGAGATTAGAATAGTGG
AGGAGGAGTGTGGTCTTGGGAACTTAGTACTCAGGAAGGAGAGGAAAAAACAAGCCATTAGGGAAAA
ATACAACCTGAAGATAGGAAACCTACGTGAAGCAGCAGAGAAAAACATATTCACCTGGATTGTTGAAGAA
GTTAAGAAAGTTCTTTGAAAATATGGCATTCTTTTCT

ACGCGTACGCGGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)


Protein Sequence: >RR214039 representing NM_001024328
 Red=Cloning site Green=Tags(s)

MAELDDDSLRLVLVGKTGSGKSATANTILGQKIFTSRIAPHAVTQTCQKASRRWKERDLLVVDTPGLFDT
 KVNLETTSEISRCVLQSCPGPHAILVLQLNRFTIEEQETVTRIKAFGKAVMKYLIILFTRKDELEDQ
 NLNDFIEDSDTNLKSIIKECDSDRYLAINNKAEGAEGEMQVQELMGFVESLVRNGLYFSDPIYKYAEQR
 LKKQVGILREIYTDVLEKEIRIVEEECGLGKLSLQEGEEKIQAIKEKYNLKNLREAAEKNIFTWIVEE
 VKKVLLKIWHFFS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001024328

ORF Size: 879 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001024328.2, NP_001019499.1

RefSeq Size: 1297 bp


RefSeq ORF: 882 bp

Locus ID: 500113

Cytogenetics: 4q24

MW: 33.4 kDa

Product images:



Circular map for RR214039