

## Product datasheet for **RR214053L3V**

### Stx3 (NM\_031124) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Stx3 (NM_031124) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Stx3
Synonyms:	Stx3a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031124
ORF Size:	867 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214053).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031124.1</a> , <a href="#">NP_112386.1</a>
RefSeq Size:	1054 bp
RefSeq ORF:	870 bp
Locus ID:	81802
UniProt ID:	<a href="#">Q08849</a>
Cytogenetics:	1q43
Gene Summary:	nervous-system specific protein; may be a component of the SNARE complex; involved in docking and fusion of synaptic vesicles to the plasma membrane; crucial for exocytosis in the nervous system [RGD, Feb 2006]



[View online »](#)