

Product datasheet for **RR214062L3V**

Gch1 (NM_024356) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gch1 (NM_024356) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gch1
Synonyms:	Gch
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024356
ORF Size:	723 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214062).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024356.1 , NP_077332.1
RefSeq Size:	1016 bp
RefSeq ORF:	726 bp
Locus ID:	29244
UniProt ID:	P22288
Cytogenetics:	15p14
Gene Summary:	an enzyme for tetrahydrobiopterin (BH(4)) synthesis; may mediate cell death of striatal neuronal and non-neuronal cells [RGD, Feb 2006]



[View online »](#)