

## Product datasheet for **RR214080L3V**

### Comp (NM\_012834) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Comp (NM_012834) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Comp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012834
ORF Size:	2265 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214080).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012834.1</a> , <a href="#">NP_036966.1</a>
RefSeq Size:	2421 bp
RefSeq ORF:	2268 bp
Locus ID:	25304
UniProt ID:	<a href="#">P35444</a>
Cytogenetics:	16p14
Gene Summary:	extracellular matrix protein that interacts with fibronectin and may play a role in extracellular matrix assembly and molecular interactions [RGD, Feb 2006]



[View online »](#)