

Product datasheet for **RR214082**

Fam83a (NM_001130576) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fam83a (NM_001130576) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fam83a
Synonyms:	RGD1565606
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR214082 representing NM_001130576 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCCGGTCAAGGCATGTGGCAAGATCCGGAACGTCTGGAAGATGTCAAGAACCAGTGGATCCGGC
CAGCCAGAGCTGATTTTAGTGATAATGAGAGCGCCCGCTGGCCACAGATGCCCTCTGGATGGGGTCC
AGAGGCTTATTGGCGGCTCTCAGCCAGGAAGGAGAGGTAGACTTCTGTCTCGGCAGAGGCCAGTAC
ATTCAGGATCAGGCCAAGGAGCCCCGGTGCCCCAGATCTTGCAGGAGGACCGAGTCAGGACCCAGAG
GACTGGACTCCTGCTCTTACAATCAGGTACCTACTTCCCTGTGGCTTCGGAAGGCAGTGAACCTGCTCT
GCTACACACCTGGACCTTGGCTGAGAAGCCCTACATGAAGGAAAAGTCCACCGCCACTGTCTACTCCAG
ATGGACAAACACAACAACATCAGAGACTTGGTCCGCCGCTGCATCTCCCGGGCCAGCCAGGTCTGGCCA
TCCTGATGGATGTGTTACGGATGTGGAGATCTTTGTGACATCTAGAGGCAGCCAACAAGCGTGGGGT
ATTTGTCTGTGTCTTCTGGATCAGGGAGGTGTGAAGCTTCCAGGAGATGTGTGACAAAATCCAGATC
TCTGATATTCATCTCAAGAACATCTCCATCAGGAGTGTGGAGGGGAGGTGTACTGTGCCAAGTCAGGCA
GGAAATTTTCGGCCAAATCCAAGAGAAGTTCATCATCTCAGACTGGAGTTTGTCTTATCCGGATCATA
CAGCTTCTCTGGCTCTGTGGACAGTGCACAGGAACATCCTCTCCAAGTTCACGGGCCAGGCGGTGGAG
CTGTTTGTGAGGAATTCGACGCCTCTACGCCTCCTCCAAGCCCTGATGGGTCTGAAGTCCCCGAGGC
TGGTCGCCCCACTCCAATCCAGCAAGGGCCAGAAAGCCCCAATGGCCGCTCAGTGGCAGCAGTGACTC
TGCCAGTGACCGCACCTCCTCAACCCCTCAGCAGCCTGTCGACTGGCAGCAACACCCACAACCGGAGT
CTATCCACGTCTCAGGGCCGGCAGTCCCCTGGCCCCATCCCACCGTGGTTCGAAGGCTGCAGCCCT
ATCACAGTACCCGGAGAGCCGTGGCCCCACAGGATCTCTTGGTCCCATGAGACCCATAAGGGAGCACC
AGCTCCCTACAGCAGCTCATGGCTACAGGCCACCCGGCTGCAGCTGCAGGAGCTCGGCTGGTGCC
AGGGTACTCAACCTGGAGCCTTTTCTACAGGCCTTACCTCATTTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR214082 representing NM_001130576
 Red=Cloning site Green=Tags(s)


MSRSRHVGKIRKRLLEDVKNQWIRPARADFSDNESARLATDALLDGGPEAYWRALSQEGEVDFLSSAEAQY
 IQDQAKEPPGAPDLAGGTESGPRGLDSCSLQSGTYFPVASEGSEPALLHTWTLAEKPYMKEKSTATVYFQ
 MDKHNNIRDLVRRICISRASQVLAALMDVFTDVEIFCDILEAANKRGVFCVLLDQGGVKLFQEMCDKIQI
 SDIHLKNIIRSVEGEVYCAKSGRKFSGQIQEKFIIISDWRFLVSGYSFVLCGHVHRNLSKFTGQAVE
 LFDEEFRRLYASSKPLMGLKSPRLVAPLQSSKGPAPNGRLSGSSDSASDRTSSNPFSSLSTGSNTHNRS
 LSTSSGPGSPLAPIPPVPRLLQPYHSTRRAVAPQDLLVPMRPHKGAPAPYSSLMAYRPTLQLQELGLVP
 RVTQPWRPFLQALPHF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001130576

ORF Size:	1308 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001130576.1 , NP_001124048.1
RefSeq Size:	1473 bp
RefSeq ORF:	1311 bp
Locus ID:	500873
Cytogenetics:	7q33
MW:	48.3 kDa