

Product datasheet for **RR214102L4V**

Npy5r (NM_012869) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Npy5r (NM_012869) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Npy5r
Synonyms:	NPYR5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012869
ORF Size:	1335 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214102).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012869.1 , NP_037001.1
RefSeq Size:	2475 bp
RefSeq ORF:	1338 bp
Locus ID:	25340
UniProt ID:	Q63634
Cytogenetics:	16p14
Gene Summary:	receptor for neuropeptide Y and neuropeptide Y Y1; may be involved in regulation of energy homeostasis [RGD, Feb 2006]



[View online »](#)