

## Product datasheet for **RR214103L3V**

### Adam17 (NM\_020306) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Adam17 (NM_020306) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Adam17
Synonyms:	TACE
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020306
ORF Size:	2481 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214103).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020306.2</a> , <a href="#">NP_064702.1</a>
RefSeq Size:	4126 bp
RefSeq ORF:	2484 bp
Locus ID:	57027
UniProt ID:	<a href="#">Q9Z1K9</a>
Cytogenetics:	6q16
Gene Summary:	expression occurs in response to oxygen glucose deprivation, may play a role in inhibition of apoptosis via induction of TNF-alpha release [RGD, Feb 2006]



[View online »](#)