

## Product datasheet for **RR214138L3V**

### Cd6 (NM\_175577) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cd6 (NM_175577) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cd6
Synonyms:	OX52
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_175577
ORF Size:	1995 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214138).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_175577.3</a> , <a href="#">NP_783167.2</a>
RefSeq Size:	2947 bp
RefSeq ORF:	1998 bp
Locus ID:	25752
Cytogenetics:	1q43
Gene Summary:	cell surface molecule found on T lymphocytes that may regulate CD5 tyrosine phosphorylation; associates with CD5, Lck, Fyn, ZAP-70 and Itk [RGD, Feb 2006]



[View online »](#)