

Product datasheet for **RR214144**

Asb17 (NM_001109676) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Asb17 (NM_001109676) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Asb17
Synonyms:	asb-17
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR214144 representing NM_001109676 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAATAACTCTTCTAAATTATGCCGGAAGACTTCTTTTCCAAGAGGCAATATATTCTGCAACCTCATTG
ACAAGATAGTTAAGAGGCCCTCTTGCAGTTTCTGGGTCAGTGGGGTACCACTGCTATGAACCCAGGAT
TTACAGAACGCTGGCAAAAATTCTGAGGTATGTTGACTTGGATGGCTTTGACATACTCCTCACTGACTAT
ATTGCTTTTGTGGAAAAGTCGGGACACCGTTTCGAACTCAATTTAACCTCGAGTTTACTGAAATATGTG
TGAATACCATTCTGTATTGGGTTTTTGGCAGAAAAGGTAACCCTGACTTTGTGGAAGTCTCTCAAGAA
GACGAAGGACTATGTCCAAGACAGAAGCTGCAGCCTAGCACTCATATGGAGAACCCTCACACCCGTGTAC
TGCCCCAGTCCCCTGAGTGGCATCACACCTCTACTCTATGTGGCTCAGACAAGACAGTCAAATATCTTAA
AAATTCTCCTGCAGTATGGAATCCTAGAAAAGAGAAAAAACCTATCAACATTGTCCTGACAATACTACT
TTACCCCTCGAGAGTGAGAATAATGGTTGATCACGAGTTGATTGACATTCAGAAGATGCCAAGACGTGT
TAAACGTTATGTTCCAGAGTGCTTTCTGTAATCTCAGTCAGGGAGATAGAGACACAGCTGAACCTTAGGAA
GACGCCCAATTATCAAATTGGTTGGACTACATCCCTCCAACAAGATACAAAGATCCATGTGAACCTCGT
CCACCTCTGCAGAATAACCATTAGAACCAACTGCTGGCCAACAATATGCTCCCAATGGAATATTTTCG
CTTCTAATTCCTACTCGTCTACAAAACCTCCTGAATTTAGAAAAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR214144 representing NM_001109676
 Red=Cloning site Green=Tags(s)

MNNSKLCRKTSFPRGNIFCNLIDKIVKRPSLQFLGQWGYHCYEPRIYRTLAKILRYVDLDGFDILLTDY
 IAFVEKSGHRFELNFNLEFTEICVNTILYWVFARKGNPDFVELLLKTKDYVQDRCSLALIWRFTPVY
 CPSPLSGITPLLVAQTRQSNILKILLQYGILEREKNPINIVLTILLYPSRVRIMVDHELIDIQEDAKTC
 LTLCSRVL SVISVREIETQLNLGRRP IIQNWLDYIPTRYKDPCELVHLCRITIRTQLLANMLPNGIFS
 LLIPTRLQNFLNLES


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001109676

ORF Size: 885 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001109676.2 , NP_001103146.1
RefSeq Size:	930 bp
RefSeq ORF:	888 bp
Locus ID:	687364
Cytogenetics:	2q45
MW:	34.5 kDa