

Product datasheet for RR214150

Prap1 (NM_031669) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Prap1 (NM_031669) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Prap1
Synonyms: Upa
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR214150 representing NM_031669
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGAGGTTTCTCCTGGCCACCTGTCTGGTGGCTGTGCTGCTGTGGGAGGCAGGTGCAATTCAGCAC
ACCAGGTCCTGTCAAGACCAAAGGCAAACATGTGTTCCCTGAACAGGAAACAGAGAAGGCCTGGGGTAC
CAGAGCCATGGAACCTTTGAAAAGGATGACCAGCTGAGGGCACTGCTTCTGTACCCAAGCAGAAGCTT
GCAGCCACTGAGGAAAAGCATTAGATACCATGACCTGGGTGGAACTAAGGACATCCTGAGCCGTTTTTC
GGAATCCTCTTACAGGTCCAGAAGCTGGATCTTGACAGCATATACCATCCCATGTCTGAGGATGTCCAGAA
TGAAGAAGTGCCCGAGTCTAGGCCATCTATATCGTCAGGTGCTGCACGGACCAGAGGAGGACCTTGAC
CACATCTCCCACTCTCTGGAGGATTCTGGAGAGCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR214150 representing NM_031669
Red=Cloning site Green=Tags(s)

MKRFLLATCLVAVLLWEAGAIPAHQVPVKTKGKHVFPEQETEKAWGTRAMEPLEKDDQLRALLPVPKQKL
AATEEKHSDMTWVETKDILSRFRNPLQPELDLDSIYHPMSDEVQNEEVPQSRPILYRQVLHGPEEDLD
HISHSLEDSGEP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_031669

ORF Size: 456 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_031669.1</u> , <u>NP_113857.1</u>
RefSeq Size:	671 bp
RefSeq ORF:	459 bp
Locus ID:	60574
UniProt ID:	<u>Q9ES75</u>
Cytogenetics:	1q41
MW:	17.3 kDa
Gene Summary:	specifically expressed in intestine and uterus; may play a role in pregnant uterus, and in digestion and absorption of nutrients in the gastrointestinal tract [RGD, Feb 2006]