

Product datasheet for **RR214150L4V**

Prap1 (NM_031669) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Prap1 (NM_031669) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Prap1
Synonyms:	Upa
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031669
ORF Size:	456 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214150).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031669.1 , NP_113857.1
RefSeq Size:	671 bp
RefSeq ORF:	459 bp
Locus ID:	60574
UniProt ID:	Q9ES75
Cytogenetics:	1q41
Gene Summary:	specifically expressed in intestine and uterus; may play a role in pregnant uterus, and in digestion and absorption of nutrients in the gastrointestinal tract [RGD, Feb 2006]



[View online »](#)