

Product datasheet for **RR214151L3V**

Ace2 (NM_001012006) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ace2 (NM_001012006) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ace2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001012006
ORF Size:	2415 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214151).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001012006.2 , NP_001012006.1
RefSeq Size:	3268 bp
RefSeq ORF:	2418 bp
Locus ID:	302668
UniProt ID:	Q5EGZ1
Cytogenetics:	Xq14
Gene Summary:	enzyme peptidyl-dipeptidase; mouse homolog has a role in regulation of cardiac contractility and regulating Angiotensin 1 and 2 levels [RGD, Feb 2006]



[View online »](#)