

Product datasheet for **RR214158L4V**

Ednra (NM_012550) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ednra (NM_012550) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ednra
Synonyms:	Endor; ET-A; ET-AR; Eta; RGD1559432
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012550
ORF Size:	1278 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214158).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012550.2 , NP_036682.1
RefSeq Size:	4595 bp
RefSeq ORF:	1281 bp
Locus ID:	24326
UniProt ID:	P26684
Cytogenetics:	19q11
Gene Summary:	binds endothelin 1 and induces an increase in intracellular calcium [RGD, Feb 2006]



[View online »](#)