

## Product datasheet for **RR214192L4V**

### Fem1a (NM\_001025706) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Fem1a (NM_001025706) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fem1a
Synonyms:	MGC114562
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001025706
ORF Size:	1962 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214192).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001025706.1</a> , <a href="#">NP_001020877.1</a>
RefSeq Size:	2222 bp
RefSeq ORF:	1965 bp
Locus ID:	316131
UniProt ID:	<a href="#">Q4V890</a>
Cytogenetics:	9q12
Gene Summary:	Probable component of an E3 ubiquitin-protein ligase complex, in which it may act as a substrate recognition subunit. May participate in antiinflammatory signaling via its interaction with PTGER4 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)