

Product datasheet for **RR214193L4V**

Lenep (NM_053614) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lenep (NM_053614) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Lenep
Synonyms:	Elp; Lep503
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053614
ORF Size:	183 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214193).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053614.1 , NP_446066.1
RefSeq Size:	593 bp
RefSeq ORF:	186 bp
Locus ID:	113917
UniProt ID:	Q9WVB7
Cytogenetics:	2q34
Gene Summary:	lens epithelial cell gene; may be involved in lens epithelial cell differentiation [RGD, Feb 2006]



[View online »](#)