

Product datasheet for **RR214204L4V**

Tcf12 (NM_013176) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tcf12 (NM_013176) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tcf12
Synonyms:	SCBPA; SCBP alpha
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013176
ORF Size:	2121 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214204).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013176.2 , NP_037308.1
RefSeq Size:	2531 bp
RefSeq ORF:	2124 bp
Locus ID:	25720
UniProt ID:	P51514
Cytogenetics:	8q24
Gene Summary:	basic helix-loop-helix (bHLH) E-protein that recognizes the consensus binding site CANNTG (E-box); may be a transcriptional activator [RGD, Feb 2006]



[View online »](#)