

Product datasheet for **RR214209L4V**

Arl6ip1 (NM_198737) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arl6ip1 (NM_198737) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Arl6ip1
Synonyms:	MGC116042
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_198737
ORF Size:	609 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214209).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198737.2 , NP_942032.1
RefSeq Size:	2139 bp
RefSeq ORF:	612 bp
Locus ID:	293551
Cytogenetics:	1q35
Gene Summary:	mouse homolog interacts with ADP-ribosylation-like factor-6 (ARL6); may play a role in hematopoietic maturation [RGD, Feb 2006]



[View online »](#)