

Product datasheet for **RR214217L3V**

Nedd4 (NM_012986) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nedd4 (NM_012986) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nedd4
Synonyms:	Nedd4a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012986
ORF Size:	2667 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214217).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012986.1 , NP_037118.1
RefSeq Size:	5373 bp
RefSeq ORF:	2670 bp
Locus ID:	25489
UniProt ID:	Q62940
Cytogenetics:	8q24
Gene Summary:	may regulate function of the amiloride-sensitive epithelial sodium channel ENaC [RGD, Feb 2006]



[View online »](#)