

Product datasheet for **RR214267L3V**

Kiss1r (NM_023992) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kiss1r (NM_023992) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Kiss1r
Synonyms:	Gpr54
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_023992
ORF Size:	1188 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214267).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023992.2 , NP_076482.2
RefSeq Size:	1502 bp
RefSeq ORF:	1191 bp
Locus ID:	78976
UniProt ID:	Q924U1
Cytogenetics:	7q11
Gene Summary:	a G protein-coupled receptor; human homolog may bind ligand and metastasis suppressor gene KiSS-1 [RGD, Feb 2006]



[View online »](#)