

Product datasheet for **RR214314L3V**

Pxn (NM_001012147) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pxn (NM_001012147) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pxn
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001012147
ORF Size:	1758 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214314).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001012147.1 , NP_001012147.1
RefSeq Size:	3737 bp
RefSeq ORF:	1761 bp
Locus ID:	360820
UniProt ID:	Q66H76
Cytogenetics:	12q16
Gene Summary:	Cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion).[UniProtKB/Swiss-Prot Function]



[View online »](#)