

## Product datasheet for RR214355L3

### PPP2R3A (NM\_001012202) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PPP2R3A (NM_001012202) Rat Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	PPP2R3A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214355).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



RefSeq ORF: 3453 bp  
Locus ID: 363122  
Cytogenetics: 8q31