

## Product datasheet for **RR214362**

### Nsun7 (NM\_001017452) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Nsun7 (NM_001017452) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nsun7
Synonyms:	RGD1310137
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR214362 representing NM_001017452 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGTTGGATCCGAACAGTGCAATGGACTTTTTGACCAAGAAAATATGGAGGAGCTCTCCAACTCGCTT  
CCCTAGAGATGTCAGGGGATTTAGTTTCCAAGAGAAAACCTTCAGCTGCGCTAGAAAATCTAGCTATCC  
AGACTCCGTGTACGTTACGGCGGCAACATTTTTCAGGGTATTCGAATCCAAAGGCCAGCAGACAAAGTC  
ATCATAAATTACGGGGATGAACTCTTGCCGTCCTCTGGACCCAGGCCGAGGAAGATTCTTTTCAGCGTT  
TGTCTATGAGTTGGCTTTCAGTACCCTGAAGTATCAAGATATTTTGAAAGTATATTAATAGATACCGA  
TATCTTCTCAAGTACCACAATCCAGATCAGTTGAACAGCCTTATCATTGTGATGCTGTATGATTTCCAA  
GATAGAAAATTTAGCCCCGTGTCCTTCCCAATGACGCAGTGACCATACCAGAAGTTCAAGAAGTAGAGA  
ACCTTCTTAACAGCTTTAAGACAAAGTTGGCTGCGGCATTGGCGAGATGTCGGATCAAGCATGATGCCCT  
TTCGATTTACCATATTCTGCCAGAGACAGTTAGGAAACAGGAACAGCGGGCCTCGACTTTGCCTCTGTAT  
GCTTGGATAAATATTGGTAAAGTAAGCCTTGAAGAAGTTACGATGATCTGAGGAAAAAGGTTTTGATA  
AGGTCAAGTCGATAGCGTATCTCAATGACAAGGCTTTTGCCGTTGATCAGCATTGCTTCAATGTCTTAAT  
TTTTCCATCTCATTTCAAAGAGAACTTTTAAACATAGATCTTATAAAAGACTATAAACTCATTTTCCAG  
GACAAATCTCGAAGCCTTGCTGTCCACTGTCAAGGCTTAAATAAACATGGACGACGATATCCTAATGG  
TCAACACAGGCTCATGGTACACAGTTGCCCATATGTCCATCTTAACAAGTGGCAACACCTCGAAAATATT  
TGTATGTGGAATAGAATCGGAAGAGAAGGATTCTGATGCGAAGAACTTTTCACAAAAATGGGATGTGAA  
AGTATGTGCCGACGTCTATTTGACCAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR214362 representing NM\_001017452  
 Red=Cloning site Green=Tags(s)


MLDPNSAMDFDQENMEELSQLASLEMSGDLVSKRKPSAALEKSSYPDSVYVTAANIFQGIRIQRPADKV  
 IINYGDELLPSSWTQAEEDSFQRLSYELAFSTLKYQDILESILIDTDFSTTIPDQLNSLIIVMLYDFQ  
 DRKFQPRVLPNDVAVTIPEVQEEVENLLNSFKTKLAAALARCRIKHDALSIYHILPETVRKQEQRASTLPLY  
 AWINIGKVSLEEVYDDLKKGFDKVKSIAYLNDKVFVAVDQHCFNVLIFPSHFKRELLNIDLKDYKLIFQ  
 DKSRSLAVHSVKALINMDDDILMVNTGSWYTVAHMSILTSGNTSKIFVCGIESEEKSDAKKLFCKMGCE  
 SMCRRFLFDH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001017452

<b>ORF Size:</b>	1077 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001017452.1</a> , <a href="#">NP_001017452.1</a>
<b>RefSeq Size:</b>	1436 bp
<b>RefSeq ORF:</b>	1080 bp
<b>Locus ID:</b>	305339
<b>Cytogenetics:</b>	14p11
<b>MW:</b>	41 kDa