

Product datasheet for **RR214380L3V**

Prickle1 (NM_199396) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Prickle1 (NM_199396) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Prickle1
Synonyms:	Rilp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_199396
ORF Size:	2493 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214380).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_199396.1 , NP_955428.1
RefSeq Size:	4485 bp
RefSeq ORF:	2496 bp
Locus ID:	315259
Cytogenetics:	7q35
Gene Summary:	plays a role in the nuclear localization of transcriptional regulators REST4 and REST/NRSF; may act as a nuclear receptor for transcription factors [RGD, Feb 2006]



[View online »](#)