

Product datasheet for **RR214386**

Synpr (NM_023974) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Synpr (NM_023974) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Synpr
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR214386 representing NM_023974
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGTGTATGGTGATATTTGCTCCGCTTTTTGCAATCTTTCGATTGCAACATGTGGCGGCTATCCGGAG
GCCTGCGGCTGAGCGTGGACTGCGTCAACAAGACAGAAAGTAACCTCAGCATAGACATAGCCTTTGCATA
TCCATTCAGGCTGCACCGGTGACCTTTGAAGTGCCACGTGTGAGGAAAGGAACGGCAGAAGCTGGCG
CTGGTTGGAGACTCCTCATCCTCAGCAGAATTCTTTGTCACCGTGGCTGTGTTTCGCTTTCTACTCCC
TGGCTGCCACTGTCGTGTATATTTCTCCAGAACAAGTACCGCGAAAACAACCGGGGCCACTCATTGA
CTTTATCGTCACTGTAGTCTTTTCATTCTTGTGGCTGGTGGGTTTCGCTGCCTGGGCAAAAGGACTGTCT
GATGTCAAAGTGGCAACAGATCCCAAAGAAGTCTGTTGCTGATGTGAGCTTGAAGCAGCCTTCCAACA
AGTGCATGGCTGTCCATAGCCCTGTCATGTCCAGCTTAAATACATCTGTGGTCTTCGGGTTTTGAACTT
TATCCTCTGGGCTGGAACATTTGGTTTTGTTTTCAAGGAGACTGGCTGGCATTCTCGGGACAGAGATAC
CTTTCCGACCCCATGGAGAAGCACTCCAGCAGCTACAATCAAGGCGGCTACAACCAAGACAGTTATGGGT
CAAGTGGTGGGTACAGTCAGCAGGCGAGTCTGGGCAACTTCTGATGAGTTTGCCAACAGCCTAGTGG
TCCCACTTCTTTAACAATCAAATA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_023974.1 , NP_076464.1
RefSeq Size:	1586 bp
RefSeq ORF:	798 bp
Locus ID:	66030
UniProt ID:	P22831
Cytogenetics:	15p16-p15
MW:	29.1 kDa
Gene Summary:	synaptic vesicle membrane protein; member of the synaptophysin family [RGD, Feb 2006]