

## Product datasheet for **RR214391**

### **Akap8l (NM\_001013946) Rat Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Akap8l (NM_001013946) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Akap8l
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RR214391 representing NM\_001013946  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCCGCATCGCC

ATGAGCTACACAGGCTTTGTCCAAGGATCCGAAACCCTTTGCAGTCTACCTACTGTGACACCAGTGCCC  
AGCCACCTGTGATTATGGATATGGAACCTTGGAACTCTGGGACAAATAGAGGCTACGAGAACTATGGTTA  
TGGCTATGGCTATGGCCAGGATAACACCACCAACTATGGGTATGGTATGGCCACTTCACACTCTTGGGAA  
ATGGCTAACTCTGACACAAATGCAAACCCTAGTGCCTCGGGTAGCGCCAGTGCCGATTCCGTTTTATCCA  
GAATTAACCAGCGCTTAGATATGATGCCACACTTGAGACAGACATGATAACAAGGAGGTGTGTACGGCTC  
AGGTGGTGAGAAAGATATGATTCCTATGAGGCCTGTGACTCAAGGGCCATCCTGAGTGAACGAGACCTC  
TACCGGTCGAGCTATGACTATGGCGAGCTTGACCCTGAGATGGAAATGGCCTATGAGGGTCAATATGATG  
CCTACCGTGACCAGTTTGAATGCGTGGAGGAGATACATTTGGTCCAAGGGCTCAGGGCTGGGCCGGGA  
TGCCCGGAGTGGACGGCCAATGGCTCGGGCTATGGGCGCATGTGGGAAGACCCATGGGGGCCGGGGC  
CAGTGCATGCCTGGTGCCTCCCGCTGCCTCCCTCTTCTCCAGAACATCATACCCGAGTATGGCATGT  
TCCAGGGCATGCGTGGTGGGGTGCCTTCTGCTGGTGGCTCCCGATTTGGCTTCGGATTTGGCAATGGCAT  
GAAGCAGATGAGGCGGACCTGGAAGACCTGGACCACAGCGGACTTCCGGACAAAAGAAGAAGAGAAAAG  
CAGGGTGGCAGTCTGATGAGCCAGACTGCAAGGCTACCAGAACAGACTGCTCAGACAACAGTGACTCAG  
ACAACGATGAAGTACTGAGGGGGAAGCTGCAGAGGGTGTGAAAGTGTGAGGCTCTGGAGAAAGGCTC  
CCGAGCAGAAGGAGAGGACGAGGAGGGAAAAGAAGATGGGAGAGAAGAAGGCAAAGAAGATTTGGAGAAA  
GGGGCCCTGACCGCCAGGATGAGAGCAGCCAGGCCAAGCGCAAGTTGCAGGCAAGCAAGAAGAGCCAAG  
ACAAACAGAAAAAGCGGCAGCGGACCGTATGGTGGAAAGGATCCAGTTTGTCTGTCTCTGCAATA  
CCGGACCTTCTATGAAGACGAGATGGGTAGCCACCTTGATAGCAAGTTCCACAAGGAACACTTCAAATAT  
GTAGGCACCAAGCTTCCAAGCAGACAGCTGACTTTCTGCAGGACTATGTCACCAACAAGACCAAGAAGA  
CAGAGGAGCTCCGAAAACTGTAGAGGACCTTGATGGTCTGATCCAACAATCTACAGAGACCAAGATCT  
GACCCAAGAAATTGCTATGGAGCATTTTGTGAAGAAGGTGGAGGCGGCCATTGTGCAGCCTGTGACCTT  
TTCATCCCATGCAGTTTGAATCATCCAAAAGCATCTCAAGACTATGGATCACAACCGAAACCGCAGGC  
TCATGATGGAGCAGTCCAAGAAGTCCTCGCTTATGGTGGCACGCAGTATTCTGAACAACAAGCTCATTAG  
CAAGAAGCTGGAGCGATATCTTAAGGGCGAGAATCCTTTTACTGACAACCCGGAGGAGGAGAAGGAACAA  
GACGAGGTGGAGGCTGGCGCTCTGGAAGAAGGGGCACCAAGCGAAGCGACAGGGCTCACAGAAGGCGTGT  
CGGTCCAGCCCCGGTGCCTTGGAGCCACCCCGGGGACTGCAACACCGCCACCTCCGCCACCACCAGA  
GGAGGAGAACCAGTACCCTTGCTGGGCGGGGCTACAGTGCCAGATCCGCGGGATCCCGGGGCTCGAC  
GTAGAAGACGACGAGGAGGGCGGTGGGGTCTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR214391 representing NM\_001013946  
Red=Cloning site Green=Tags(s)


MSYTGfVQGSETTLQSTYCDTSAQPTCDYGYGTWNSGTNRGYENYGYGYGYGDNTTNYGYGMATSHSWE  
MANSDTNANPSASGSASADSVL SRINQRLDMMPHLETDMIQGGVYSGGGERYDSYEACDSRAILSERDL  
YRSSYDYGELDPEMEMAYEQYDAYRDQFRMRGGDTFGPRAQQWARDARSGRPMASGYGRMWEDPMGARG  
QCMPGASRLPSLFSQNI IPEYGMFQGMRRGGAFSGGSRF GFGNGMKQMRRTWKTWTTADFRTKKKKRK  
QGGSPDEPDCKATRTDCSDNSDSDNDEGTEGEAAEGAESA EALEKGSRAEGEDEEGKEDGREGKEDLEK  
GALTAQDESSQAKRKLQASKKSQDKQKKRQRDRMVERIQFVCSLCKYRTFYEDEMGSHLDSKFHKEHFKY  
VGTKLPKQTADFLQDYVTNKTKKTEELRKTVEDLDGLIQQIYRDQDLTQE IAMEHFVKKVEAAHCAACDL  
FIPMQFGI IQKHLKTM DHNRRRLMMEQSKSSLMVARSILNNKLI SKKLERYLKGENPFDNPEEEKEQ  
DEVEAGALEEGAPSEATGLTEGVSVQPPVPLEPTPGTATPPPPPPPEEENPVPLLGALQCQIRGIPGLD  
VEDDEEGGGP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**


**ACCN:** NM\_001013946

**ORF Size:** 1923 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001013946.1](#), [NP\\_001013968.1](#)  
RefSeq Size: 2194 bp  
RefSeq ORF: 1926 bp  
Locus ID: 299569  
Cytogenetics: 7q11  
MW: 71.3 kDa