

Product datasheet for **RR214393L3V**

Ppara (NM_013196) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ppara (NM_013196) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ppara
Synonyms:	PPAR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013196
ORF Size:	1404 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214393).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013196.1 , NP_037328.1
RefSeq Size:	2022 bp
RefSeq ORF:	1407 bp
Locus ID:	25747
UniProt ID:	P37230
Cytogenetics:	7q34
Gene Summary:	nuclear steroid receptor; heterodimerizes with retinoid-X receptor (RXR), and binds to peroxisome proliferator response elements (PPREs) to initiate the transcriptional regulation of target genes [RGD, Feb 2006]



[View online »](#)