

Product datasheet for **RR214404**

Rpl24 (NM_022515) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rpl24 (NM_022515) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rpl24
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR214404 representing NM_022515
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGGTCGAGCTGTGCAGTTTTAGTGGGTACAAGATCTACCCGGGACACGGGCGGCGCTACGCCAGGA
CCGATGGGAAGGTTTTCCAGTTTTCTTAATGCCAAATGTGAGTCCGCATTCCTTTCCAAAAGAACCCCTCG
GCAAATTAACCTGGACTGCCTCTACAGAAGAAAACACAAGAAAGGACAGTCGGAAGAAATTCAAAAGAAA
AGAACCAGCCGTCAGTCAAGTTCCAGCGGCCATCACAGGCGTCTCTGGCTGATATAATGGCCAAGA
GGAATCAGAAACCAGAAGTTAGGAAAGCTCAGCGAGAACAGGCTATCAGGGCTGCCAAGGAAGCAAAAAA
GGCTAAGCAGGCATCAAAGAAGACAGCAATGGCTGCTGCCAAGGCTCCCACAAAGGCAGCCCTAAACAA
AAGATTGTGAAGCCTGTGAAGTCTCTGCTCCAGAGTTGGTGGGAAACGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR214404 representing NM_022515
Red=Cloning site **Green**=Tags(s)

MKVELCSFSGYKIYPGHGRRYARTDGKVFQFLNAKCESAFLSKRNPRQINWTVLYRRKHKKGQSEEIQKK
RTRRAVKFQRAITGASLADIMAKRNQKPEVRKAQREQAIRAAKEAKKAKQASKKTAMAAKAPTKAAPKQ
KIVKPKVSAPRVGGKR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN:

NM_022515

ORF Size:

471 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_022515.1</u> , <u>NP_071960.1</u>
RefSeq Size:	541 bp
RefSeq ORF:	474 bp
Locus ID:	64307
UniProt ID:	<u>P83732</u>
Cytogenetics:	11q12
MW:	17.8 kDa
Gene Summary:	60S ribosomal subunit protein [RGD, Feb 2006]