

## Product datasheet for **RR214424**

### **Pdp1 (NM\_019372) Rat Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Pdp1 (NM_019372) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pdp1
Synonyms:	Ppm2c
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)



**Cloning Scheme:**



**ACCN:** NM\_019372

**ORF Size:** 1614 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_019372.4](#), [NP\\_062245.2](#)

**RefSeq Size:** 2515 bp

**RefSeq ORF:** 1617 bp


**Locus ID:** 54705

**Cytogenetics:** 5q13

MW: 61.2 kDa

Gene Summary: catalytic subunit of the mitochondrial matrix phosphatase; involved in reactivation of the pyruvate dehydrogenase complex [RGD, Feb 2006]

### Product images:



Circular map for RR214424