

Product datasheet for **RR214444L3V**

Sparcl1 (NM_012946) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sparcl1 (NM_012946) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sparcl1
Synonyms:	Ecm2; Sc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012946
ORF Size:	1902 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214444).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012946.1 , NP_037078.1
RefSeq Size:	2676 bp
RefSeq ORF:	1905 bp
Locus ID:	25434
UniProt ID:	P24054
Cytogenetics:	14p22
Gene Summary:	putative extracellular matrix glycoprotein that binds calcium and plays an important role in the developing brain [RGD, Feb 2006]



[View online »](#)