

Product datasheet for **RR214462L3V**

Zfp143 (NM_001012169) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Zfp143 (NM_001012169) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Zfp143
Synonyms:	STAF; Znf143
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001012169
ORF Size:	1878 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214462).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001012169.1 , NP_001012169.1
RefSeq Size:	3040 bp
RefSeq ORF:	1881 bp
Locus ID:	361627
UniProt ID:	Q5XIU2
Cytogenetics:	1q33



[View online »](#)

Gene Summary:

Transcriptional activator. Activates the gene for selenocysteine tRNA (tRNA^{Sec}). Binds to the SPH motif of small nuclear RNA (snRNA) gene promoters. Participates in efficient U6 RNA polymerase III transcription via its interaction with CHD8 (By similarity).[UniProtKB/Swiss-Prot Function]