

## Product datasheet for **RR214465L3V**

### **Ccn1 (NM\_053662) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ccn1 (NM_053662) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ccn1
Synonyms:	Ania-6a; Ccn1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053662
ORF Size:	1581 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214465).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053662.3</a> , <a href="#">NP_446114.1</a>
RefSeq Size:	2128 bp
RefSeq ORF:	1584 bp
Locus ID:	114121
UniProt ID:	<a href="#">Q9R1Q2</a>
Cytogenetics:	2q31
Gene Summary:	cyclin protein that is regulated by cAMP response element binding protein (CREB); may be involved with differential RNA processing, such as alternative splicing [RGD, Feb 2006]



[View online »](#)