

Product datasheet for **RR214478**

Rnase4 (NM_020082) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rnase4 (NM_020082) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rnase4
Synonyms: Rab1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR214478 representing NM_020082
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACATACAGAGGCCAATCCTTGCTCTGCTCTTGTGCTGACCCTGCTGGGTTAGGGCTTGTAC
AGCCCTCCTATGGCCAAGATAGAATGTACCAACGGTTCCTTAGACAGCATGTGGACCCTGAGGGGACAGG
CGGCAGCGACAACACTACTGCAACGTGATGATGCAGAGACGGAGGATGACTTCTACCCAGTGCAAACGCTTC
AACACCTTCATCCACGAAGACATCTGGAACATTCGCAGCATCTGTGATACTGCCAATATCCCATGCAAGA
ATGGCAATATGAACTGTACGAAGGCATAGTGAGGGTCACTGACTGCAGAGAGACAGGGAGCTCTGTGCC
CCACAACGTAGGTACAGGGCGAGAGCCAGCACTAGGCGAGTTGTATTGCCTGTGAGGGTACCCAGAG
GTCCCAGTGCACTTTGACAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR214478 representing NM_020082
 Red=Cloning site Green=Tags(s)

MDIQRQTSLLLLLLLTLLGLGLVQPSYGQDRMYQRFLRQHVDPEGTGSDNYCNVMMQRRRMTSTQCKRF
NTFIHEDIWNIRSICDTANIPCKNGNMNCHEGIVRVTDCRETGSSVPHNCRYRARASTRRVVIACEGTP
VPVHFDR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_020082

ORF Size: 441 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_020082.2</u> , <u>NP_064467.1</u>
RefSeq Size:	1360 bp
RefSeq ORF:	444 bp
Locus ID:	56759
UniProt ID:	<u>O55004</u>
Cytogenetics:	15p14
MW:	16.9 kDa
Gene Summary:	member of a family of pyrimidine-specific ribonucleases that displays a substrate preference for poly(U) over poly(C) [RGD, Feb 2006]