

Product datasheet for **RR214519**

TP53RK (NM_001108606) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: TP53RK (NM_001108606) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: TP53RK
Synonyms: Trp53rk
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR214519 representing NM_001108606
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGGCGCTCTTCGGAGGCGGAGGCGGAGGCGCTGGCCGCGGCGCGAGCGGAGCCGCTTCC
TGAGCGGCCTGGAGCTGGTGCAGCAGGGCGCCGAGGCTCGAGTGTCCGTGGCCGCTTCCAGGGCCGAGC
GGCCGTGGTGAAGCACCGCTTCCCTAAGGGCTACCGGCACCCGGAGCTGGAGGCGCGCTCGGCCGCGG
CGGACGGTTCAGGAGGCTCGCGCGCTGTGCGCTGCCGCCGTGCGGGAATAGCTGCCCGCTCGTCTTCT
TTGTGGACTATGCTTCTAACTGCCTATATATGGAAGAAATCGAAGACTCGGTGACTGTTGAGATTATAT
TCAATCCACTATGGAGACTAAACAAGACCCCAATGCCTCTTGGACTTGGCCAGGAAGATCGGGCAGGTT
CTGGCCAGACTGCATGACGAGGACCTCATTACGGGGACCTCACCACCTCCAACATGCTCCTGACACGGC
CCCTGGAGCAGCTGCACATCGTGCTCATCGACTTCGGGCTGAGTTTTGTCTCAGGACTGCCAGAAGACAA
AGGCGTAGACCTCTATGCTCCTGGAGAAAGCCTTCTCAGCACGCACCCCACTGAGACTGTGTTTGA
GCCTTTCTGAAGATTACGGTGCCTCATCAAGAAGTCCGATCCAGTGTGAAGAAGTTAGATGAGGTGC
GCCTGAGAGGGCGAAAGCGGTCCATGGTTGGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR214519 representing NM_001108606
 Red=Cloning site Green=Tags(s)

MAGASSEAEAEALAAARERSRLFLSGLELVQQGAEARVFRGRFQGRAAVVKHRFPKGYRHPELEARLGRR
 RTVQEARALLRCRRAGIAAPVVFFVDYASNCLYMEEIEDSVTVRDYIQSTMETKQDPQCLLDLARKIGQV
 LARLHDEDLIHGDLTTSNMLLTRPLEQLHIVLIDFGLSFVSGLPEDKGVLDLVLEKAFSLTHPHTETVFE
 AFLKSYGASSKSDPVLKKLDEVRLRGRKRSMVG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001108606

ORF Size: 732 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001108606.1](#), [NP_001102076.1](#)

RefSeq Size: 2312 bp


RefSeq ORF: 735 bp

Locus ID: 362272

Cytogenetics: 3q42

MW: 27.4 kDa

Product images:



Circular map for RR214519