

## Product datasheet for **RR214547L4V**

### **Ndufs6 (NM\_019223) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ndufs6 (NM_019223) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ndufs6
Synonyms:	lp13dis; Ndub13; RATlp13dis
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019223
ORF Size:	348 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214547).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019223.1</a> , <a href="#">NP_062096.1</a>
RefSeq Size:	351 bp
RefSeq ORF:	351 bp
Locus ID:	29478
UniProt ID:	<a href="#">P52504</a>
Cytogenetics:	2q12
Gene Summary:	subunit of NADH:ubiquinone oxidoreductase [RGD, Feb 2006]



[View online »](#)