

Product datasheet for **RR214549L3V**

Rbm45 (NM_153306) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rbm45 (NM_153306) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rbm45
Synonyms:	Drb1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153306
ORF Size:	1428 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214549).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153306.1 , NP_695218.1
RefSeq Size:	1778 bp
RefSeq ORF:	1431 bp
Locus ID:	266631
UniProt ID:	Q8CFD1
Cytogenetics:	3q24
Gene Summary:	RNA-binding protein with binding specificity for poly(C). May play an important role in neural development.[UniProtKB/Swiss-Prot Function]


[View online »](#)