

Product datasheet for **RR214555**

Sebox (NM_023951) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sebox (NM_023951) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Sebox
Synonyms:	Og9x
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR214555 representing NM_023951 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGCCAGTCCTGTGGAGGCATCTCCAGGCTGTGCCAGTGGGTTGGTCCCCATCGGAGGAAGAGGACCA
 CTTTGTAGTTGGGCAGTTGCTGGAGCTGGAACGGGTATTTGCAGCTAGGCCCTATCCTGACATCAGCAC
 CCGTGAGCACCTGGCTCAGATAACTCACCTGCCAGAAGCAAGATCCAGGTGTGGTTCCAGAACCGTCGA
 GCCAAGAGAATCAAGGACAGAAAGCCAGGAGCCCTAAACGCCAGGCTGGAGCTTCCCCAACTGTTCTC
 TTCCAGACACTCCCAGCTATCTTGGGACCCTGGAACACCAAGCCATCCTCTACACTCTACCAGTTCAGC
 TCAGCAGACTTCGGCATGCCACACAGACCTCCTGCCTAGCCCCATTTTGGGTCCAGGGCAGAGTTGG
 TCAGGGGCTAAAGCTGCAGCCCCATGGGGGACAAGTGAGGTGTCAGGAGTTCACCTCTTAGAGCAAAATTG
 TTCTCAGACTTCACTGGGCAACCTGTCTGACCTATCTATACCTCAGCCATCGTCACCAATGTAGACCA
 CTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RR214555 representing NM_023951 Red=Cloning site Green=Tags(s)
-------------------	--

MASPVEASPGCASGLGPHRRKRTTFSVGQLELERVFAARPYPDISTREHLAQITHLPEAKIQVWFQNR
 AKRIKDRKPGALNARLELPPNCSLPDTPQLSWDPGTPSHPLHSTSSAQQTSAQPPQTSCLAPILGPGQSW
 SGAKAAAPWGTSEVSGVHSLEQIVPQTSGLNLSLDLIYTSIAIVTNVDHS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_023951

ORF Size: 564 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_023951.1](#), [NP_076441.1](#)

RefSeq Size: 567 bp


RefSeq ORF: 567 bp

Locus ID: 58922

UniProt ID: [Q9ERS8](#)

Cytogenetics: 10q25
MW: 20.2 kDa
Gene Summary: paired-like homeobox containing gene; mouse homolog is maternally expressed [RGD, Feb 2006]

Product images:



Circular map for RR214555