

Product datasheet for **RR214563L3V**

Nr5a2 (NM_021742) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nr5a2 (NM_021742) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nr5a2
Synonyms:	Ftf; Lrh-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_021742
ORF Size:	1680 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214563).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021742.1 , NP_068510.1
RefSeq Size:	1931 bp
RefSeq ORF:	1683 bp
Locus ID:	60349
UniProt ID:	Q9QWM1
Cytogenetics:	13q13
Gene Summary:	an orphan nuclear receptor that may bind DNA and activate gene transcription [RGD, Feb 2006]



[View online »](#)