

## Product datasheet for **RR214566L3V**

### Ilf3 (NM\_053412) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ilf3 (NM_053412) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ilf3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053412
ORF Size:	2730 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214566).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053412.1</a> , <a href="#">NP_445864.1</a>
RefSeq Size:	3168 bp
RefSeq ORF:	2733 bp
Locus ID:	84472
UniProt ID:	<a href="#">Q9JIL3</a>
Cytogenetics:	8q13
Gene Summary:	nuclear protein that interacts with and regulates the activity of protein-arginine methyltransferase I [RGD, Feb 2006]



[View online »](#)