

Product datasheet for **RR214581**

Msl1 (NM_001107048) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Msl1 (NM_001107048) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Msl1
Synonyms:	RGD1564778
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR214581 representing NM_001107048
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGACCATGAGATCCGCAGTGTTCAGGCGGCCGCGGCCCTGCCGGCGGCAACCCGAGCAGCAGTGG
 ACTACGAGCGGGCTGCGGCTCTGGGCGGGCCGAGGACGAGTCCGGGGCGGCCGAAACCCACTTCCTCCC
 CCGGCATCGTAAGCTCAAGGAGCCCGGGCCCCGCTGGCCTCTACCCAGGGCGGGAGCCCCGCGCCTCT
 CCAGCCGGCTGCGGCGCGGCAAGGGCCGGGTTTGTACTCCCGGCGGGGCGGCCCGGGCAGCAGG
 AAGAGAGCTGGGCGGTTTCGGTGCCTTGCCCTGTCCGCCCCGGCTACCAACAAGCCGCATCGGCGG
 GGAGCCAGTCGCAGCCGGCGCTGGCTGCAGCCCCCGGCCAAGTATCAGGCGGTGCTGCCATTAGACG
 GGCTCTCTCGTGGCGGCGCCAAAGAGCCTACGCCCTGGGCTGGGACAAGGGTGGGCGGCTCCCCAG
 CTGCCACCGCCTCGGACCCGGCGGGACCCCACTACTCTGCCCGGGCGGCCACCCCTCGCGCCAC
 CGCCACTGCCGGGACCTGGCGGCCAGTGAGGCGAGATGGAAGAGTATAAGGAAGAGCCCTCTCGGGGT
 GCGGGCGGCTCGGAGCCTCCAGTCAGGCGCCTGCCTCAAACAGATCCTTCTGCTGCAATTGGACCTCA
 TCGAACAGCAGCAGCAGCAGTTCGAAGCCAAGGAAAAGGAGATCGAGGAGCTGAAGTCGGAGAGAGATAC
 GCTCCTTGCTCGGATTGAACGTATGAAAAGGCGGATGCAGCTGGTAAAGAGGGATAACGAGAAGGAAAGG
 CACAAGCTGTTTCAGGGCTATGAACCTGAAGAGAGAGAGGAAACGGAGTTGTCTGAGAAAATTAAGTGG
 AGCGCCAGCCGGAGCTCTGCGAGACATCCAGGCTCTGCCTTCTAAGCCTTTCATGTGGCCGGAGTGG
 CAAGGGACACAAAAGGAAACCCATTTGGAATACAGAAAGAAAGACTCCTGTTAAAAAGCTGGCTCCT
 GAATTTTCAAAATCAAAACAAAACCTCTAAGCACTCTCCATTAAAGAGGAGCCCTGTGGTTCATAT
 CAGAACTGTTTGTAAACGTGAATTGAGGAGCCAAGAAACCCAGAAAAGCCCGGTCTTCAGTGGATAC
 CCCACCAAGACTCTCGACTCCCCAAAAGGGACCCAGCACCCACCCCAAGGAGCAAGCCTTCTCAAGTGAG
 ATGGAAGATTTGCCGTACCTTTCCACCACAGAAATGATTTGTGTCGTTGGCACCAGCCTCCCCATCAC
 CGTTACCATTACGGGAATCCTCTCAAAGAAGGAGGAGACTGTAGCAATTCCTTCTTGAGGGACCACTC
 AGTAGAGCCTCTAAGGGACCCAAATCCTTCAGACATTTTGAGAACCTGGATGACAGTGTGTTTTCAAAG
 AGGCATGCGAACTGGAGCTAGACGAGAAGAGGAGAAAACGATGGGATATTCAGAGGATCAGGGAACAAA
 GAATTTTACAGCGACTGCAGCTCAGAATGTATAAAAAGAAAGGAATTCAGGAATCTGAGCCTGAGGTTAC
 CTCATTTTCCCTGAGCCCGATGACGTTGAAAGTTTGATGATTACCCCTTCTTGCTGTGTAGCATTT
 GGACGGCCATTACCAAAATTAACCTCCACAGAACTTTGAACTACCTGGTTGGATGAGCGAAGCCGTTGCA
 GATTGGAGATCCAGAAGAAGCAACACCTACCGGACTGTAGGAAG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR214581 representing NM_001107048
 Red=Cloning site Green=Tags(s)


MTMRSVFKAAAAAPAGGNPEQRLDYERAAALGGPEDESGAAETHFLPRHRKLKEPGPPLASTQGGSPAPS
 PAGCGGKGRGLLLPAGAAPGQQEESWGGSVPLPCPPATKQAGIGGEPVAAGAGCSRPKYQAVLPQT
 GSLVAAAKEPTPWAGDKGAAPPAATASDPAGPPPLPLPGPPPLAPTATAGTLAASEGRWKSIRKSPLGG
 GGGSGASSQAACLKQILLQLDLIEQQQQQLAKEKEIEELKSERDTLLARIERMERMQLVKRDNEKER
 HKLFQGYEPEEREETELSEKIKLERQPELCETSQALPSKPFSCGRSGKGHKRTPFGNTERKTPVKKLAP
 EFSKFKTKTPKHSPIKEEPCGSISETVCKRELRSQETPEKPRSSVDTPPRLSTPQKGPSTHPKEQAFSSE
 MEDLPYLSTTEMYLCRWHPPPSPLPLRESSPKKEETVAIPSWRDHSVEPLRDPNPDIENLDDSVFSK
 RHAKLELDEKRRKRWDIQRIREQRILQRLQLRMYKKKGIQSEPEVTSFFPEPDDVESLMTFPLPVVAF
 GRPLPKLTPQNFELPWLDERSRCLIEIQKQTPHRTCRK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001107048

ORF Size: 1797 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001107048.1](#), [NP_001100518.1](#)

RefSeq Size: 2768 bp


RefSeq ORF: 1800 bp

Locus ID: 303514

Cytogenetics: 10q31

MW: 65.8 kDa

Product images:



Circular map for RR214581