

Product datasheet for **RR214582L3V**

Lax1 (NM_001017491) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lax1 (NM_001017491) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Lax1
Synonyms:	MGC108827
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001017491
ORF Size:	1218 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR214582).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001017491.1 , NP_001017491.1
RefSeq Size:	2241 bp
RefSeq ORF:	1221 bp
Locus ID:	498232
UniProt ID:	Q5FVQ5
Cytogenetics:	13q13
Gene Summary:	Negatively regulates TCR (T-cell antigen receptor)-mediated signaling in T-cells and BCR (B-cell antigen receptor)-mediated signaling in B-cells.[UniProtKB/Swiss-Prot Function]


[View online »](#)